mediline®

PRODUCT FROM THE GCE GROUP

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- Health care equipment – oxygen therapy and emergency devices*
- regulators, flow-meters and accessories for medicinal gases
 - humidifiers, breathing masks and canals
 - suction
 - transport ventilators and bag valve masks
 - oxygen concentrators

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*For further information see www.gce.cz.
Health Care Equipment catalogue available on request.*

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*Products and their features are in accordance with standards of Czech Republic and Slovak Republic only.
Other versions are available on request – please ask us for appropriate article numbers.*

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Medidave pressure regulators are designed to reduce the pressure of medical gas from cylinders. They are diaphragm single stage regulators with an inlet pressure gauge and are available in several versions. The regulators are fitted with an inlet filter to guard against impurities, and have a safety valve to stop the rise in working pressure. Inlet connections are presented in the table below in accordance with domestic standards (Czech and Slovak republics). Other versions are available on customer's request. Medidave regulators are designed to fulfil requirements of EN 738-1.

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The 100 regulator has a gauge to show the inlet pressure, and has a preset outlet pressure with an adjustable outlet flow. The regulator body is made of chrome plated brass. The inlet connection is fitted with a hand tightened nut. The outlet is fitted with a float-type flowmeter with a tolerance of ±10%. The regulator is designed as a direct oxygen source for patients. In long-term oxygen therapy, it is recommended to humidify oxygen (e.g. with Mediwet 200 humidifier).

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The 300 regulator has a gauge to show the inlet pressure, and has a preset outlet pressure. The regulator body is made of chrome plated brass. The inlet connection is fitted with a hand tightened nut. The outlet connection is fitted with a medical quick release coupler, in accordance with standards of Czech and Slovak republics. The regulator is designed as a supply source for medical equipment (anaesthetic machines, ventilators, ambulance pipelines, etc.).

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The 500 regulator has a gauge to show the inlet pressure, and an adjustable outlet pressure, and gauge to show the outlet pressure or flow. The regulator body is made of chrome plated brass. The inlet connection is fitted with a hand tightened nut. The regulator is designed as a versatile source either direct to a patient or for other medical equipment. In long-term oxygen therapy, it is recommended to humidify oxygen (e.g. with Mediwet 200 humidifier).

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The pressure regulator has a gauge to show the inlet pressure, and an adjustable outlet pressure, and a gauge to show the outlet pressure. The regulator body is made of brass, powder painted (black coloured), inlet and outlet connection pieces are made of chrome plated brass. The inlet connection is fitted with a hand tightened nut. The outlet connection piece is fitted with a pressure relief valve with a manual control device to easy release of the gas from a connected equipment, and with a ball closing valve. The Fixicontrol+ CO₂ regulator is suitable e.g. for filling of bags for dry gas carbon baths.

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The Mediline Demand valve resuscitator provides a fast and simple and effective mean to ventilate non-breathing patients and provides 100% oxygen to patients suffering from respiratory distress. Mediline demand valve resuscitator comes standard with an audible pressure relief system. The conveniently placed sensitive manual ventilation button on the front of the demand valve helps to ensure a patient airway by enabling both hands to remain on the mask to facilitate superior mask to face seal. The manual flow rate is fixed at 40 l/min. During spontaneous breathing the Mediline Demand valve resuscitator is designed to provide 100% oxygen to a breathing patient with minimal respiratory effort. The oxygen demanded flow is controlled by the patient's inspiratory efforts with flow rates up to 160 l/min. The Mediline Demand valve incorporates 360° patient valve swivel housing and is easily disassembled for cleaning.

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The Mediline demand valve provides a safe and effective mean of delivering 100% oxygen by inhalation, in line with the flow rate demanded by the patient. Lightweight, portable and extremely durable the Mediline demand valve is designed for the demands of the pre-hospital and in-hospital environment and can be used anywhere a regulated oxygen supply is available. The Mediline demand valve incorporates 360° patient valve swivel housing and is easily disassembled for cleaning.

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The Mediline demand valve provides a safe and effective mean of delivering 100% oxygen by inhalation, in line with the flow rate demanded by the patient. Lightweight, portable and extremely durable the Mediline demand valve is designed for the demands of the pre-hospital and in-hospital environment and can be used anywhere a regulated oxygen supply is available. A range of Mediline demand valves with adapters is offered to fit the most commonly available Bag Valve Mask resuscitators allowing an easy way to provide 100% oxygen.

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One mask fit all sizes from infants through to children. For infants and children the mask is rotated 180°. The inner-seal lip edge allows the mask to offer superior mask to face seal even with a beard or moustache. There are also children, infant and neonatal masks available.

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The scavenger adapter is easily fitted to the Mediline Demand valve to enable venting of gases when so is wished.

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Mediflow and Mediselect pressure regulators are designed to reduce the pressure of medical gasses from cylinders. They are diaphragm single stage regulators with an inlet pressure gauge, and an adjustable outlet flow. They are fitted with an inlet filter to protect against impurities, with a safety valve to stop excessive rise in the working pressure. Inlet connections presented in the table below are in accordance with the domestic standards (Czech and Slovak Republics). Other versions are available on customer's request. Mediflow and Mediselect regulators are designed to fulfil requirements of EN 738-1.

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The Mediflow regulator has an inlet pressure gauge, a preset outlet pressure and an adjustable outlet flow. The regulator body is made of nickel plated brass. The inlet connection has a nut that can be hand tightened. The outlet is fitted with a float-type flowmeter with a precision $\pm 10\%$. The regulator is designed as a direct oxygen source for patients. In case of long term oxygen therapy, it is recommended to humidify oxygen (e.g. with a Mediwet 250 humidifier).

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The Mediselect regulator has an inlet pressure gauge, a preset outlet pressure and an adjustable outlet flow. The regulator body is made of nickel plated brass. The inlet connection has a nut that can be hand tightened. The outlet is fitted with a flowmeter with a precision $\pm 10\%$. The adjustments can be made with an operating hand wheel in 10 steps:

Mediselect 1 (MC4101): 0 - 0,1 - 0,2 - 0,3 - 0,4 - 0,5 - 0,6 - 0,7 - 0,8 - 1 l/min

Mediselect 5 (MC4105): 0 - 0,5 - 0,75 - 1 - 1,5 - 2 - 2,5 - 3 - 4 - 5 l/min

Mediselect 15 (MC4115): 0 - 1 - 2 - 3 - 4 - 5 - 6 - 9 - 12 - 15 l/min.

The regulator is destined as a direct oxygen source for patients. In a case of a long term oxygen therapy, it is recommended to humidify oxygen (e.g. with a humidifying bottle Mediwet 250).

► Redukční ventily řady Mediselect, Medicontrol / *Mediselect, Medicontrol Pressure Regulators*

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Mediselect and Medicontrol pressure regulators are designed to reduce the pressure of medical gasses from cylinders. They are diaphragm single stage regulators, with an inlet pressure gauge. Mediselect regulators have an adjustable outlet flow, Medicontrol regulators have a preset outlet pressure and outlet flow. Both regulators are fitted with an inlet filter to protect against impurities, and a safety valve to stop excessive rise in the working pressure. Inlet connections presented in the table below are in accordance with domestic standards (Czech and Slovak Republics). Other versions are available on customer's request. Mediselect and Medicontrol regulators are designed to fulfil requirements of EN 738-1.

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Mediselect regulators have an inlet pressure gauge, a preset outlet pressure, and an adjustable outlet flow. The regulator body is made of nickel plated brass. The inlet connection has a nut that can be hand tightened. The outlet is fitted with a flowmeter with a tolerance of ±10%. The adjustments can be made with hand operating wheel in 10 steps: 0 - 1 - 2 - 3 - 4 - 5 - 6 - 9 - 12 - 15 l/min. The Hospistar is designed as a direct oxygen source for patients. In case of a long term oxygen therapy, it is recommended to humidify oxygen (e.g. with a Mediwet 250 humidifier). The Mediselect 15S is additionally equipped with a quick release coupling outlet (preset pressure and flow) to supply other equipment or central piping (e.g. in an Ambulance).

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The Medicontrol has a gauge to show inlet pressure, and a preset outlet pressure and flow. The regulator body is made of nickel plated brass. The inlet connection has a nut that can be hand tightened. The outlet is fitted with a quick release coupler in accordance with domestic standards (Czech and Slovak Republics). The regulator is designed as a supply source for medical equipment (anaesthetic machines, ventilators, ambulance piping etc.).

► Průtokoměry řady Mediflow, Debistar, Medimeter *Mediflow, Debistar, Medimeter flow-meters*

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Mediflow, Debistar and Medimeter gas flow-meters are design to regulate and measure the flow of medicinal gases. They are equipped with an inlet filter to protect against impurities. The inlet connection is fitted with a quick probe made of stainless steel, which enables to connect the flow-meter to a medicinal gas source (pressure regulator, wall outlet,...). The Mediflow, Debistar and Medimeter flow-meters are designed to fulfil requirements of EN 13220, EN 738-1, EN 12218, EN 739.

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Mediflow 15 (Debistar 15): 0 - 1 - 2 - 3 - 4 - 5 - 6 - 9 - 12 - 15 l/min.

Mediflow 50 (Debistar 15): 0 - 5 - 10 - 15 - 20 - 25 - 30 - 35 - 40 - 50 l/min.

Mediflow (Debistar) flow-meter regulates the flow of medicinal gases in ten preset steps using an operating hand wheel. The adjusted flow value is readable in a window on the edge of the operating hand wheel. The flow-meter's body is made of nickel plated brass. The inlet connection is fitted with a quick probe for connecting to a source of medicinal gas (pressure regulator, wall outlet). Outlet connection can be connected to a breathing mask, nasal canulas or humidifying bottle in a case of long term oxygen therapy (e. g. Mediwet 250).

The Mediflow (Debistar) flow-meter is available in oxygen and air version and with flow values: 0 - 1 l/min, 0 - 5 l/min, 0 - 15 l/min, 0 - 50 l/min.

Mediflow 1 (Debistar 15): 0 - 0,1 - 0,2 - 0,3 - 0,4 - 0,5 - 0,6 - 0,7 - 0,8 - 1 l/min.

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The Medimeter flow-meters are float-type flow-meters. They are designed to regulate the flow of medicinal gases using an operating hand wheel. The adjusted flow value is readable in a vertical position of the float. The flow-meter's body is made of chrome plated brass, upper plastic parts are made of autoclavable polycarbonate. The inlet connection is fitted with a quick probe for connecting to a source of medicinal gas (pressure regulator, wall outlet). Outlet connection can be connected to a breathing mask, nasal canulas or humidifying bottle in a case of long term oxygen therapy (e. g. Mediwet 200, Mediwet 250).

The Medimeter flow-meter is available in oxygen and air version and with flow values 0,1 - 1 l/min, 0,4 - 4 l/min, 0,5 - 5 l/min, 1,2 - 12 l/min, 1,5 - 15 l/min, 3 - 30 l/min.

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Mediwet humidifying bottles are designed to increase a relative humidity of medical oxygen supplied to a patient in a case of long term oxygen therapy. The humidifiers can be fitted to oxygen pressure regulators, oxygen flow meters, oxygen concentrators or other medicinal oxygen sources. They are filled with sterile water.

Plastic parts of the humidifier are made of polycarbonate resistant to a mechanical stress, connectors are made of chrome plated brass. The humidifiers are autoclavable at 134°C

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Breathing masks and nasal canulas are designed for supplying a self-breathing patient with medicinal oxygen. They are made of plastic, disposable. The connection hose has a special internal safety profile preventing from the oxygen delivery interruption in a case of twisted or pressed hose.

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Medievac vacuum regulators are designed to regulate vacuum in suction systems connected to a central vacuum piping. The regulators are available in 3 versions: maximum vacuum value -250 mbar, -600 mbar, -1000 mbar. The regulator is equipped with a vacuum-meter, a control hand wheel to adjust the vacuum value and an ON-OFF lever for immediate opening-closing of the vacuum flow. The regulator can be fitted to a wall outlet directly or can be fixed on a rail and be fitted to the wall outlet using a flexible hose.

The regulator's body is made of shock-proof plastic, the connections are made of chrome plated brass. The inlet quick probe is made of stainless steel.

It is recommended to use the regulator fitted with a safety suction jar Medicollect 300 (it contains an over-load protection valve to protect the regulator and the central vacuum piping from a contamination); it is recommended to connect the regulator to collecting suction jars Medicollect using a silicone suction hose.

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The MS-32 venturi ejector is an alternative source for suction where no central vacuum piping is available. The ejector creates vacuum by connecting to a source of regulated pressure air or pressure oxygen (wall outlets, pressure regulators). The ejector is fitted with a vacuum-meter, a control hand wheel to adjust the vacuum value and an ON-OFF button for immediate opening-closing of the vacuum flow. The ejector also contains a noise filter for the appropriate operator comfort. It is recommended to use the ejector with Medicollect collecting suction jars; it is recommended to connect them using a silicone suction hose.

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The Mediair venturi ejector is an alternative source for suction where no central vacuum piping is available. The ejector creates vacuum by connecting to a source of regulated pressure air (wall outlets, pressure regulators). The ejector is fitted with a vacuum-meter and a control hand wheel to adjust the vacuum value. The ejector also contains a noise filter for the appropriate operator comfort. It is recommended to use the regulator fitted with a safety suction jar Medicollect 300 (it contains an over-load protection valve to protect the regulator and the central vacuum piping from a contamination); it is recommended to connect the regulator to collecting suction jars Medicollect using a silicone suction hose.

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The Mediclic venturi ejector is a simple alternative source for suction where no central vacuum piping is available. The ejector creates vacuum by connecting to a source of regulated pressure air or oxygen (wall outlets, pressure regulators). The ejector is fitted with a control hand wheel to adjust the vacuum value.

It is recommended to use the ejector with Medicollect collecting suction jars; it is recommended to connect them using a silicone suction hose.

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Medicollect 300 safety suction jar is suitable for use with Medievac vacuum regulators or Mediair venturi ejectors. It is fitted with an overflow valve to protect the vacuum source and the vacuum central piping against a contamination. Medicollect safety jar is autoclavable at 134°C. It is available with or without a bacterial filter. The filter is a subject to replace each use.

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Medicollect collecting suction jars are suitable for use in hospital environment. They should be connected to a vacuum source using silicone suction hoses. They are fitted with an overflow valve to protect the vacuum source and the vacuum central piping against a contamination. The jars are made of polycarbonate with chrome plated brass nipples. They are autoclavable at 134°C.

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Disposable suction bags Medibag are designed for collecting of liquid materials in hospital environment. They are disposable – must be destroyed after filling in accordance with hospital regulations. When using, the Medibag is fixed in a special designed jar with a lid.

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